

FlexiFuge

BioTool Microcentrifuge

Put a new Spin on Things

FlexiFuge is a versatile, personal micro-centrifuge for quick low speed centrifugation. Whisper quiet motor offers two speeds (2000g or 5000g / 6,000 or 9,200 RPM) depending on application. FlexiFuge can be operated in continuous mode, pulsed for a short burst of rapid centrifugation, or timed to spin from 1-15 minutes. Regardless of setting, rotors stop once the lid release button is pushed.



A ring of LED light illuminates when the unit is active and can be viewed through the transparent lid. Colour can be changed (red, green, blue) or turned off by a switch under the unit. LED Lights glow Red, Green or Blue or can be turned off by switch on bottom of unit.



FlexiFuge is supplied with an 8-Place fixed angle rotor for 15mL tubes, a rotor for (4), 8 x 0.2mL PCR tube strips and 8 adapters for 0.5mL tubes. It can also accommodate a 4-place rotor (purchased separately) to spin 5mL micro-centrifuge tubes,



- AG.T2076, clear
- AG.T2077, blue
- AG.T2078, yellow
- AG.T2079, green
- AG.T2080, rose



AG.C1000	FlexiFuge, complete with 3 Rotors, 220/240V power supply	1 ea
AG.C0300	Spare Rotor for 1.5ml Microcentrifuge Tubes (1 inclusive)	1 ea
AG.C0301	Spare Rotor for 0,2ml PCR Tube-Strips (1 inclusive)	1 ea
AG.C0302	Special Rotor for Argos Exclusive 5ml Snap-Cap Tubes	1 ea
AG.C0303	Spare Rotor for 0.5ml Microcentrifuge Tubes (1 inclusive)	1 ea
AG.T207□	5ml Snap-Cap Tubes, PP, (withstands - 90°C up to 120°C/15min, 15psi)	250 / bag