

# Why You Need the MiniSpace® and Select IUI® Catheters

Both the MiniSpace® and the Select IUI® insemination catheters are designed specifically for your patients. The Luer-lock hub allows convenient syringe attachment.

## The Select IUI

The dual openings of the Select IUI mirror the normal orientation of the Fallopian tubes, while the wings of the Luer-lock hub indicate the position of the openings. With an effective length of 17cm, the Select IUI has an interior volume of 0.34cc and an outer diameter of 2.1mm for the distal 4.5cm of length.

## The MiniSpace

The ultra-low residual volume MiniSpace leaves next to nothing in the catheter. With a single opening at the tip and an effective length of 17.2 cm, the MiniSpace has an interior volume of 0.07cc and an outer diameter of 1.6mm for the distal 4.5cm of length.

The carefully designed Select IUI and MiniSpace catheters feature the finest craftsmanship but are amazingly cost effective.

