

The advantages of cryotherapy

- no destruction of biological functions and structures
- excellent cosmetic effects
- almost no scar formations
- non bleeding treatment
- self-anesthetizing
- practically painless, therefore also suitable for children

The patented "Liquid freezing" -system by N₂O

Our method is used in almost common superficial skin indications as warts, keratoses, papillomas, basilomas, precanceroses. Our cryosurgical system is as effective as therapy with N₂ and especially N₂O must not be underestimated simply because it does not reach temperatures as low as those reached by N₂. The range, not the effect within the effective range, is smaller.

Our devices provide the full freezing capacity of nitrous oxide (N₂O) from ready available mini cartridges –each cartridge with its own filter-(produced only for Cryoswiss by Linde AG). The evaporating temperature of liquid phase of -89°C is converted into therapeutic energy with almost no loss. Nitrous oxide in liquid phase is applied directly on the skin through different mini glass capillaries. After just one cryo-treatment the majority of symptoms disappear. Precise application with a constant penetrating temperature of -89°C/-128.2 F avoids application to the surrounding healthy tissue.

-89 °C
-128.2 F

... the new standard for cryosurgery!

**Liquid freezing technologies by N₂O -
the finest, most innovative simplified
cryosurgical system**

Our cryosurgical devices are CE0344 marked.
div. German patents
Chinese patent ZE 0013458.7 –
Hong Kong HK 10386886
US pat.pend.Appl.Ser.No.11/570.585,
Canada pat.pend.Appl.Ser.No.2.571.006,
Japan pat.pend.App.No.2007-518714

DISTRIBUTOR:

INTERNATIONAL SALES CONTACT:

CRYOSWISS GmbH
Rütlistrasse 57
CH-4051 BASEL
Switzerland

Answ.mach./Fax/ +41 61 271.5991
cryoalfa@cryoswiss.com



CRYOALFA PRO

- Precise dosage and exceptional simple handling controlled with user-friendly button valve
- The energy resource N₂O is sterile-filtered to 5 μ in the cartridge
- Mini-device, light & portable carrying in the overall pocket
- Various applicators are available for special applications in the gynaecology, urology, dentistry, proctology, aesthetic medicine, podiatry to be connected to PRO

CRYOALFA PERFECT

- ⇨ All devices come with 16g N₂O cartridges and 1,00 mm diameter applicator, packaged in stable case with noble skin-coloured foam
- ⇨ Accurate application to 1 mm with lever-valve
- ⇨ Suitable for a doctor's consulting room, hospital out-patient, practitioner for naturopathy
- ⇨ Several lesions or patients can be treated with one cartridge
- ⇨ Closed applicator in different sizes for contact freezing using for aesthetic dermatology can be connected

CRYOALFA LUX

- All LUX devices are working with the newest cartridge technology.
- The 5μ filtered N₂O cartridges have additionally an inbuilt valve. Now the handling is easier yet.
- Different applicators like PRO are available for special applications to be connected to Cryoalfa LUX