



THE EZ-RANGE

- EZ-Tip®
- EZ-Grip®
- EZ-Strip®
- EZ-Squeeze™

The EZ-Range of denudation and transfer devices offer unparalleled quality, from the ergonomic EZ-Grip® pipettor to the unique EZ-Strip® denudation system. RI presents an EZ solution to your professional requirements, whilst also satisfying your personal preferences. We have an option for every need.

— **High Quality**

RI EZ pipettes are rigorously checked for quality - ensuring consistent geometries, close tolerances and smooth edges.

— **Sterile and Biocompatible**

All RI EZ denudation and transfer pipettes are batch tested to ensure embryo safety using the industry’s highest pass mark in Mouse Embryo Assay (MEA) - at least 80% of the test embryos must reach the blastocyst stage. The batches are also certified free of endotoxins by the Limulus Amoebocyte Lysate (LAL) test (<0.5 Eu/ml pass mark). Our pipettes are manufactured in an ISO 13485 facility and sterilised by an approved gamma irradiation facility meeting ISO 11137 standards.

— **Sealed**

All RI denudation and transfer pipettes are available individually packaged for the best sterility guarantee. EZ-Tip® pipettes are also available in a multipack option for convenience and economy.

— **Sizes**

The EZ-Range is available in a large selection of colour coded sizes:

EZ-TIP Ø0.9

PRODUCT:	EZ-STRIP Ø1.5	MULTIPACK	SINGLE	EZ-SQUEEZE Ø1.5
CODE:	7-72-1###/20	7-72-2###/5	7-72-4###/20	7-72-3###/20
QUANTITY:	20 devices	100 pipettes	20 pipettes	20 pipettes
SIZE:	125	✓	✓	
	131	✓		✓
	135	✓	✓	✓
	145	✓	✓	✓
	155	✓	✓	✓
	170	✓	✓	✓
	200	✓	✓	✓
	290	✓	✓	✓
	600	✓	✓	✓

Take the code above (e.g. 7-72-1###/20) and add the pipette size to generate the ordering code. For example 7-72-3131/20 is the code for an EZ-Squeeze™ pack of 20, size 131µm.



EZ-Tip®

The EZ-Tip® is a premium quality plastic pipette for use with the ergonomic EZ-Grip® pipettor. EZ-Tip® pipettes are supplied either in a pouch packed vial of 20 pipettes or as individually sealed pipettes in a box of 20.



EZ-Tip® - Vial of 20



EZ-Grip®

This premium quality adjustable aspiration volume pipettor has been developed for the most comfortable and precise sample manipulation. This streamlined design is intended for practical sterilisation and long-term durability. For use with the EZ-Tip® premium quality plastic pipette.



EZ-Tip® - Singles



EZ-Grip®



EZ-Squeeze™

A combined plastic pipette and bulb system which fits to a handle accessory. The EZ-Squeeze™ comes pre-assembled - the pipette is fitted with an aspiration bulb, creating a reliable seal, which is not always guaranteed with customer assembled systems. This combined unit enables a faster set-up.



EZ-Squeeze™



Sealed EZ-Squeeze™



EZ-Strip®

EZ-Strip® is a complete single-use system - to minimise the risk of cross contamination, each EZ-Strip® is supplied sterile, with its own handle and is completely disposable.



EZ-Strip®





THE EZ-RANGE

Research Instruments Ltd, Bickland Industrial Park, Falmouth, Cornwall, TR11 4TA, UK

Tel: + 44 (0) 1326 372 753 Fax: + 44 (0) 1326 378 783 e-mail: sales@research-instruments.com www.research-instruments.com

